

The ASEPT® Pleural & Peritoneal Drainage Systems are indicated for long term, intermittent drainage of pleural effusions and malignant ascites.¹ This insertion kit comes with most items generally required for catheter implantation. The primary component of the ASEPT® Drainage System is the catheter, which has an aseptic, self-sealing valve attached that provides a secure connection. The valve can be connected to the ASEPT® Drainage Line, attached to the ASEPT® Vacuum Drainage Bottles or Drainage Bag, or a stand-alone drainage line to connect to other fluid collection equipment.

## **Features and Benefits**

## **ASEPT® Pleural & Peritoneal Drainage Systems**







ASEPT® Catheter Valve

## **ASEPT® Catheter Valve**

- SECURE CONNECTION
- **ASEPTIC DESIGN**
- NEEDLE-FREE & SWABABLE
- **EASY TO REPLACE**
- > SELF-SEALING

### **Various Drainage Options**





The product offering of the ASEPT® Drainage Kit configurations includes a 600ml or 1,000ml vacuum bottle set, a 2,000ml gravity drainage bag and the procedure pack with all the dressing change items. A stand alone drain line is available to connect to the ASEPT® Drainage Catheter for drainage with other collection equipment.

## **Accessories**







**ASEPT® Drainage Line Set** 



**ASEPT® Procedure Pack** 

Contains all the dressing change items. (Included with every ASEPT® Drainage Kit)

# **Ordering Information**

Product Description	Part Number
ASEPT® Pleural Drainage System (Includes ASEPT® drainage catheter and insertion kit)	M7001 (1 each)
ASEPT® Peritoneal Drainage System (Includes ASEPT® drainage catheter and insertion kit)	M7008 (1 each)
ASEPT® 600ml Drainage Kit (Includes vacuum bottle, drainage line and procedure pack)	M7002 (case of 10 each)
ASEPT® 1,000ml Drainage Kit (Includes vacuum bottle, drainage line and procedure pack)	M7000 (case of 10 each)
ASEPT® 2,000ml Drainage Bag Kit (Includes gravity bag with drainage line and procedure pack - for peritoneal drainage only)	M7060 (case of 10 each)
ASEPT® 2,000ml Drainage Bag - Gravity bag, for peritoneal drainage only. (Includes drainage bag with drainage line)	M7058 (case of 10 each)
ASEPT® Procedure Pack (Includes all the dressing change items)	M7040 (case of 10 each)
ASEPT® Replacement Valve	M7004 (box of 5 each)
ASEPT® Drainage Line Set (Includes drainage line with end caps, ASEPT® connector, 5-in-1 adapter and pinch clamp)	M7005 (box of 10 each)

¹ Indications for Use

The ASEPT\* Pleural Drainage System is intended for long term, intermittent drainage of symptomatic, recurrent, pleural effusions, including malignant pleural effusions and other pleural effusions that do not respond to treatment of the underlying disease. The ASEPT\* Peritoneal Drainage System is indicated for periodic drainage of recurrent and symptomatic malignant ascites. The catheter is intended for long term access of the peritoneal cavity in order to relieve symptoms such as dyspnea.

#### www.pfmmedicalusa.com



customerservice@pfmmedicalusa.com • www.pfmmedicalusa.com